

Armrests

We offer a range of different armrests. All the armrests are made of polyurethane foam and they all have a wooden core. The armrests have integrated threaded inserts which are coordinated with our armrest hardware with matched hole spacing.

Armrests, rounded



- Made of soft polyurethane flexible foam
- Easy fitting to wheelchairs and seat shells
- Easy-care and washable
- Wooden core with 100 mm hole spacing

E 1 4 100*	
A=	27,0 cm
B=	5,0 cm

Armrests, rounded, mini



- Made of soft polyurethane flexible foam
- Easy fitting to wheelchairs and seat shells
- Easy-care and washable
- Wooden core with 50 mm hole spacing

E 1 4 101*	
A=	16,5 cm
B=	5,0 cm

Armrests, flat



- Simple armrest made of polyurethane flexible foam
- Easy-care and washable
- Wooden core with 200 mm hole spacing

E 1 4 200*	
A=	30,0 cm
B=	5,0 cm

Armrests, channel form



- Channel form armrest made of polyurethane flexible foam
- Easy-care and washable
- Wooden core with 200 mm hole spacing

E 1 4 210*	
A=	30,0 cm
B=	5,0 cm

Armrests, channel-shaped, with hand support



- Anatomically shaped armrest with hand support
- This armrest made of polyurethane flexible foam caters for a wide range of treatment options
- The 5 threaded inserts are provided for easy fitting to wheelchairs and seat shells
- Available either left or right-handed
- Hole spacing in 50 mm increments

E 1 4 300R=1*	
E 1 4 300L=1*	
A=	30,0 cm
B=	11,0 cm
C=	7,0 cm

E 1 4 300R=2*	
E 1 4 300L=2*	
A=	40,0 cm
B=	13,5 cm
C=	7,0 cm

*A= Length, B= Width, C= Height

Armrests



Individual Armrest, flat, narrow

- Simple armrest made of polyurethane flexible foam
- Easy-care and washable
- Wooden core with 200 mm hole spacing


E 1 4 250*	
A=	37,0 cm
B=	7,0 cm



Individual Armrest, flat

- Simple, soft, flat armrest
- Wooden or aluminium mounting plate with chip compound foam and memory foam
- Incl. cover
- Please state measurements
- Please state measurements for mounting holes

E 1 4 401 R
E 1 4 401 L



Individual Armrest, anatomic, simple form


- Channel form straight armrest
- Wooden or aluminium mounting plate with chip compound foam and memory foam
- Incl. cover
- Please state measurements
- Please state measurements for mounting holes

E 1 4 402 R
E 1 4 402 L



- Anatomically shaped armrest for complex arm bedding/support
- Wooden or aluminium mounting plate with chip compound foam and memory foam
- Incl. cover
- Please state measurements
- Please state measurements for mounting holes

E 1 4 403 R
E 1 4 403 L



Armrest Mounting Plate

- For fitting armrests to seat shells
- Height-adjustable using two vertical slots
- Powder-coated aluminium
- Hole spacing 100 mm

E 1 4 010*	
A=	19,0 cm
B=	15,0 cm
C=	4,0 cm

*A= Height, B= Width, C= Depth