

# Assessor

The Assessor is our new versatile Assessment chair (Art. No. E 3 6 100) for the highly qualified fitting of seating shells.


This new wheelbase offers an outstanding level of adjustments and flexibility in all situations.

It offers a limitless variety of fitting options and may be used for a range of sizes from small children up to adult sizes.

**Please note the following extras are available for this article.**



**Swing away lateral supports (height adjustable)**



**Accessories Tray**



from



# ORDER FORM



McLean REHatechnik GmbH  
Breiter Anger 20  
37115 Duderstadt

Tel.: +49 (0)5527-997390  
Fax: +49 (0)5527-9973950  
E-Mail: info@mclean-rehatechnik.de





Delivery address

Com.:

## Assessor

Fitting chair for custom seating

Please note: The gas cylinder assists the tilt in space

Colour		Specification		
<input type="checkbox"/> Pearl grey	Fitting chair <ul style="list-style-type: none"><li>- ergonomic frame</li><li>- height and angle adjustment</li><li>- seat and back plates</li><li>- armrests with aluminium side plate</li><li>- back adjustment (1 gas cylinder)</li><li>- tilt in space (2 gas cylinder)</li><li>- rear whell 12" with drum brakes</li><li>- front wheels 140 mm castors</li></ul>		Seat width: 35-55 cm Seat depth: 35-60 cm Seat height: 56 cm Angle of folding back: 90° to 140° Seat plate angle: -8° to 38°	Weight: kg  Max loading 120 kg  Wheelbase width: 68 cm Total height: 118 cm
Castor wheels	Back	Seat	Accessories	
140 mm front castors Solid PU With threaded axle	Adjustable push, Handle angle, height adjustment, gas cylinder for angle, adjustment; fitted with 6 removeable back plates	Seat with aluminium baseplate and aluminium sideplates with armrest mounting. Seat angle adjustment with gas cylinder	Headrest 	Tip bar  Mounted on RH side 
Driving wheel	Foot rest	Accessories		
Driving wheel 12" with attendant operated drum brakes.	2 piece footrest with angle and height adjustment 	One piece footrest with 90° angle. (may be ordered instead of 2 piece angle adj. footrests)   This footrest is neither angle adjustable nor flip up.	Extra plastic tray in frame	Thoracic pads(flip back) with skai leather cover. With mounting plate on back