



Chest Support Harnesses

Chest Support Harnesses, made of elastic neoprene that is kind to the skin, are used with trunk instability. The elastic material offers the user greater freedom of movement. Individual adjustment of each fixing strap results in perfect adaptation to the clinical requirements. The Chest Support Harnesses are easy-care and machine-washable up to 30° C. They are available in six sizes in the standard version. Women should use the chest support harnesses designed for women. All sizes also available with a zip. These items can also be custom-made.

Größe/ Size*	A	B	C	D	E	F	Band
0 (0 R)	21,0 (21,5)	20,0 (22,0)	12,0 (14,0)	14,5 (15,5)	10,0 (11,0)	4,0 (4,0)	18 mm
0 plus (0 plus R)	27,0 (27,5)	21,0 (22,0)	13,0 (14,0)	17,5 (18,5)	11,5 (12,0)	4,0 (4,0)	20 mm
1 (1 R)	33,0 (33,5)	22,0 (23,0)	13,0 (14,0)	20,5 (21,0)	12,0 (13,0)	4,5 (4,5)	25 mm
2 (2 R)	38,0 (38,5)	25,0 (26,0)	15,0 (16,0)	23,0 (23,5)	12,5 (13,5)	5,0 (5,0)	25 mm
3 (3 R)	42,0 (42,0)	27,0 (27,5)	17,0 (17,5)	24,5 (25,0)	13,0 (14,0)	5,0 (5,0)	25 mm
4 (4 R)	48,0 (48,0)	31,0 (32,0)	20,0 (21,0)	31,0 (31,5)	16,0 (16,5)	5,5 (5,5)	25 mm
W 2 (W 2 R)	33,0 (33,0)	24,5 (26,5)	14,5 (16,5)	22,0 (23,0)	17,0 (17,5)	5,0 (5,0)	25 mm
W 3 (W 3 R)	40,0 (40,0)	26,0 (28,0)	16,0 (18,0)	25,0 (26,0)	24,0 (24,5)	5,0 (5,0)	25 mm



Chest Support Harness, mini

Art. No.	Description
E 1 5 300=0	Chest support harness, mini

Chest Support Harness



Art. No.	Description
E 1 5 300=0 Plus	Chest support harness, size 0 plus
E 1 5 300=1	Chest support harness, small
E 1 5 300=2	Chest support harness, medium
E 1 5 300=3	Chest support harness, large
E 1 5 300=4	Chest support harness, extra large

Chest Support Harness with zip



Art.-Nr.	Bezeichnung
E 1 5 300 R=0	Chest support harness, mini with zip
E 1 5 300 R=1	Chest support harness, small with zip
E 1 5 300 R=2	Chest support harness, medium with zip
E 1 5 300 R=3	Chest support harness, large with zip
E 1 5 300 R=4	Chest support harness, extra large with zip

Chest Support Harnesses

Chest Support Harness, female



Art. No.	Description
E 1 5 301 W=2	Chest support harness, female, medium
E 1 5 300 W=3	Chest support harness, female, large

Chest Support Harness, female with zip



Art. No.	Description
E 1 5 301 RW=2	Chest support harness, female, medium with zip
E 1 5 301 RW=3	Chest support harness, female, large with zip


Dynamic Torso Belt

Elastic torso belt with adjustable tension. Especially suitable for active wheelchair users with a high degree of paralysis. Developed for everyday usage and sporting activities.



Art. No.	Description
E 1 5 350=1	Dynamic Torso Belt, small
E 1 5 350=2	Dynamic Torso Belt, medium
E 1 5 350=3	Dynamic Torso Belt, large

	A	B	C
size 1	225 (200-270)	45	400
size 2	300 (270-350)	60	400
size 3	400 (370-450)	75	400



Belt Pad



Art. No.	Description
E 1 5 360	Belt Pad

- 150 x 50 mm
- unto 30 mm Belt

Chest Support Harnesses

Upper Arm Retraction Belt


A further development of the chest harness, which offer a light pull on either right, left or both shoulders.

This harness is usually used to prevent light spinal twisting or torsion which is often encountered with clients who have skoliosis.

Please state: right left both sides



Art. No.	Description
E 1 5 310=1	Upper arm retraction belt, male, size 1
E 1 5 310=2	Upper arm retraction belt, male, size 2
E 1 5 310=3	Upper arm retraction belt, male, size 3
E 1 5 310=2 W	Upper arm retraction belt, female, size 2
E 1 5 310=3 W	Upper arm retraction belt, female, size 3



	A	B	C	D	E	F
size 1	35,0	22,0	12,5	22,0	6,0	32,5
size 2	40,0	24,5	14,5	24,0	8,0	37,0
Size 3	44,0	27,0	16,0	26,0	14,0	39,0
W2	35,0	23,0	14,0	23,0	8,0	29,0
W3	42,0	24,0	14,0	26,0	11,0	33,0

*Size details are in cm. The drawing shows a female harness. Because of the great flexibility of the harnesses, the measurements given are average. Variations in size may occur.