


Pelvic Harnesses


These pelvic harnesses are used as pelvic restraints. They are made of elastic neoprene that is kind to the skin. The elastic material offers the greatest possible freedom of movement combined with secure restraint and optimum positioning of the pelvis. This soft and smooth material reduces the likelihood of pressure sores or chaffing. Individual fitting is possible by adjusting the fixing straps as required. The pelvic harnesses are easy-care and machine-washable up to 30° C. In the standard version they are available in different sizes.

Pelvic Harness





Art. No.	E 1 5 500=1	E 1 5 500=2	E 1 5 500=3	E 1 5 500=4	E 1 5 500=5	E 1 5 500=6
	size 1	size 2	size 3	size 4	size 5	size 6
A=	19,0 cm	23,0 cm	27,0 cm	31,0 cm	35,0 cm	40,0 cm
B=	35,0 cm	40,0 cm	45,0 cm	50,0 cm	55,0 cm	60,0 cm
C=	10,0 cm	12,0 cm	13,0 cm	13,0 cm	14,0 cm	14,0 cm
D=	13,0 cm	16,0 cm	19,0 cm	22,0 cm	25,0 cm	29,0 cm
E=	4,0 cm	5,0 cm	5,5 cm	5,5 cm	6,0 cm	6,0 cm

Pelvic Harness, T-Model



Art. No.	E 1 5 501=1	E 1 5 501=2	E 1 5 501=3
	size 1	size 2	size 3
A=	22,0 cm	28,0 cm	34,0 cm
B=	15,0 cm	17,0 cm	21,0 cm
C=	4,0 cm	3,5 cm	4,5 cm

Combi Pelvic Harness



Art. No.	E 15 502=1	E 15 502=2	E 15 502=3
	Size 1	Size 2	Size 3
A=	21 cm	22,5 cm	24 cm
B=	33 cm	35 cm	42 cm
C=	21 cm	29 cm	36 cm
D=	15 cm	17 cm	18 cm
E=	5 cm	5 cm	5 cm